CHIMIA-REPORT

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Leserdienst Nr. 19

The new ABI PRISM® 7900HT Sequence Detection System from Applied Biosystems

Applied Biosystems has introduced the ABI PRISM® 7900HT Sequence Detection System. This real-time PCR instrument has been

launched to speed pharmaceutical customers' drug discovery process. The 7900HT performs gene expression quantification and detection of

single nucleotide polymorphisms (SNPs) at exceptionally high-throughputs by taking a higher plate-density than previous generation instruments, without compromising speed, resolution or reproducibility.

Identifying promising targets is the first step in developing new drug compounds, traditionally a slow, tedious, and expensive process. The 7900HT Sequence Detector offers a choice of assay chemistries to speed this up. The sensitive probe-based, fluorogenic 5' nuclease assay allows multiplex reactions for real-time quantification and SNP genotyping, and the SYBR® Green 1 double-stranded DNA-binding dye assay is ideal for target identification in screening assays. Laser light is distributed to each well in the reaction plate, exciting the fluorescent dyes present. The fluorescent light generated is then directed through a spectrograph to a charge-coupled device (CCD) camera, giving spectrally resolved fluorescence data from each sample.

An exciting feature of the 7900HT is the automation accessory that loads up to eighty-four 384-well plates, for 24-hour hands-free operation and unprecedented throughput. The system also takes both 96- and 384-well plates for enhanced flexibility. Sample tracking is simplified and multiple fluorophores can be used in a single reaction.

Sample tubes remain closed throughout the PCR process to control contamination. Cost-effective, high-throughput gene expression studies are possible, as the sequence detector performs both absolute and relative quantification of target in unknown samples and quantifies both target and endogenous control in multiplex reactions. Large-scale SNP screening on the 7900HT is carried out with the fluorogenic 5' nuclease assay, with two TaqMan® probes, one for each allele, multiplexed in a single tube. The instrument then measures post-PCR endpoint fluorescence, allowing genotyping results to be available quickly.

The new ABI PRISM® 7900HT Sequence Detection System will prove invaluable for those involved in pharmaceutical research, helping to bring new drugs to market more rapidly and efficiently.

For more information on the ABI PRISM® 7900HT Sequence Detection System, please contact: abdirect@eur.appliedbiosystems.com

- APPLIED BIOSYSTEMS
Kelvin Close,
Birchwood Science Park North,
Warrington, WA3 7PB, UK
Telefon +44 (0) 1925 825650
Telefax +44 (0) 1925 282502
www.appliedbiosystems.com

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KRETZ AG